



WHAT SHOULD YOU DO WITH THE LA BASSINE POOL AFTER THE BIRTH?

1. There are disposable gloves provided for the cleaning process.
2. Remove any particulate matter and drain the tub following the directions outlined in the Instruction Guide. Use the **dark hose for removing used water**, by attaching it to the submersible pump. While draining, the pump can be placed in a colander to prevent suction to the bottom of the tub which stops removal of water. The water can drain into a shower, bathtub, or your own yard.
3. Dispose of the Disposable Liner.
4. Cleanse all surfaces of the tub and all accessories to be returned, including the floating thermometer and submersible pump, by thoroughly rinsing them with fresh water and wetting them with the cleaning solution provided (Tincture of Green Soap). Let the cleaning solution sit on the tub and any accessories for 10 minutes before wiping them down with a clean towel. Run fresh water through the pump and dark hose for at least 10 seconds.
5. Wipe all items with a dry clean towel and let air dry until completely dry. **This is an important step.**
6. Empty the hoses of any water by lifting the hose from the center and letting water drain out of the ends. Wipe down with the cleaning solution, dry, coil, and for each, screw the ends of the hose together so that no drops of water can escape.
7. Let water drain from the submersible pump, dry and place in the Red BioMed Bag provided.
8. When the tub is completely clean and dry, let the air empty, gently fold and place in the bag. When other items are completely clean and dry they can be placed in the container. Be mindful when packing the equipment to exercise care and avoid damaging the materials.

Below is a check list of all the materials to be packed for return:

- ☐ La Bassine Professional Pool & Bag (*cleaned, dried, folded, placed in bag*)
- ☐ Floating Thermometer (*drained, cleaned and dry*)
- ☐ Submersible Pump (*rinsed with fresh water, cleaned, dried, placed in Red BioMed bag*)
- ☐ Electric Air Pump and Attachments
- ☐ Spray Bottle with Remaining Cleaning Solution
- ☐ White Hose (*emptied, dried, coiled, ends screwed together*)
- ☐ Dark Hose (*rinsed with fresh water, emptied, cleaned, dried, coiled, ends screwed together*)
- ☐ Binder with Instructions
- ☐ Storage Box with Lid
- ☐ Faucet Adapter(s)
- ☐ La Bassine Patch Kit
- ☐ Any additional rented items

Call Esther at (510) 725-5456 before the return date to schedule pick up/drop off.